



Your *neffy*[®] Demonstration Device Kit

This kit contains air-filled, disposable demonstration devices designed to help familiarize you and your patients with *neffy*[®]'s administration, feel and size. The demonstration devices do not contain active drug and are intended for multiple demonstrations.

Please note: to reset this device, firmly pull the plunger back until you hear a 'click'. The device is now ready for demonstration again. Please refer to the *neffy*[®] product monograph for complete administration instructions.

We hope you'll find these demonstration devices useful when discussing the administration instructions of *neffy*[®] with patients and caregivers in your office.



SCAN HERE
to watch the *neffy*[®]
Administration Video



For further information, please visit www.neffy.ca or contact your ALK representative.

Not actual size.



ALK-Abello Pharmaceuticals Inc. All rights reserved.
neffy[®] and *neffy*[®] logo are trademarks or registered trademarks of ARS
Pharmaceutical Operations, Inc. and are used here with permission.